



CÁNULA DE GUEDEL (OROFARÍNGEA)

CÓDIGO	MEDIDA	CUADRO BÁSICO	CAJA MASTER
RSPCP656003	Natural No. 000 (30 mm)	*	50 Pzas.
RSPCP656004	Rosa No. 00 (40 mm)	060.167.3304/ 379.171.1363	50 Pzas.
RSPCP656005	Azul No. 0 (50 mm)	060.167.0458	50 Pzas.
RSPCP656006	Negra No. 1 (60 mm)	060.167.3312	50 Pzas.
RSPCP656007	Blanca No. 2 (70 mm)	060.167.0466	50 Pzas.
RSPCP656008	Verde No. 3 (80 mm)	060.167.3320	50 Pzas.
RSPCP656009	Amarilla No. 4 (90 mm)	060.167.0482	50 Pzas.
RSPCP656010	Roja No. 5 (100 mm)	060.167.3346	50 Pzas.
RSPCP656011	Naranja No. 6 (110 mm)	060.167.0680	50 Pzas.
RSPCP656012	Morada No. 7 (120 mm)	*	50 Pzas.



▪ Cánula de plástico curvo PVC FLEXIBLE con una pestaña reforzada en el extremo oral, de forma plana para asegurar que encaje perfectamente entre la lengua y el paladar duro.